PCR Cleaner™

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PCR Cleaner™

Against DNA/RNA residue

LEMNISCARE™
LemnisCare Medical Technology
Please follow the instruction
LS-SX03-1812-5001

Product Introduction

PCR Cleaner is a spray that specifically removes surface nucleic acid contaminants. It can effectively degrade nucleic acid molecules. It is easy to use and has obvious effects. It is suitable for PCR laboratories, PCR equipment and consumables.


Main Ingredients

The solution consists of a surfactant and a non-alkaline non-carcinogen, which is non-toxic and harmless to the human body. The main components include water, ethoxylates.

Function Principles

The main principle is to change the charge distribution of the nucleic acid adsorbed on the surface of the instrument and the test bench under the combined action of the effective components in the solution, to promote the detachment of the nucleic acid from the adsorbed surface, and then further combine the shear force generated during the wiping process to further contaminate the nucleic acid. Dissociated from the adsorption surface and entered the PCR Cleaner reagent system to achieve the effect of removing surface nucleic acid contamination. PCR Cleaner can non-enzymatically hydrolyze nucleic acid strands under normal temperature conditions to achieve the effect of decontamination.

Operation Instructions

Precautions

Wear safety glasses and disposable gloves when handling reagents; This product can be used in glass, ceramics, plastics, rubber, steel and metals; Avoid using PCR Cleaner on electronic devices, such as electric dispensers or balance.

FAQs

Q1: Will the reagent affect the normal experiment? A1: Does not affect the normal experiment.

Q2: What should I do if I accidentally swallow or touch my eyes? A2: After swallowing, drink plenty of water, rinse your eyes, and see a doctor in time.

Q3: I don't know if the laboratory is polluted. How to detect it? A3: You can use our nucleic acid purification kit to perform PCR detection after sampling on the refrigerator door handle, freezer frame, or test bench, centrifuge.

Q4: Is it flammable and explosive? A4: Non-inflammable and non-explosive substances, please feel free to transport and store.

Storage & Term of validity

Room temperature, avoid direct sunlight, valid for 12 months.

Product Specifications

Cat. No.	Product Name	产品规格
PCN01-0500	PCR Cleaner	Bottle A: 500 mL; Bottle B: 500 mL