Space Cleaner™

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LEMNISCARE™
LemnisCare Medical Technology
请仔细阅读说明书并按规范使用
LS-SX03-1812-5002

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LS-SX03-1812-5002

Product Introduction

Space Cleaner™ is a cleaner specially designed for pollution prevention of laboratory super clean tables. Bottle A can effectively remove DNA/RNA and other residual nucleic acid pollution, as well as eliminate false positives; Bottle B can sterilize, prevent mildew and clean the super clean table. It is suitable for super clean desk, PCR laboratory environment, PCR equipment and consumables, etc.

Main Ingredients


Bottle A: water, ethoxy compounds (surfactants) and so on; Bottle B: Ag^+ and other composite nanoionic solutions.

Function Principles


The main principle of Bottle A is to change the charge distribution of the nucleic acid adsorbed on the surface of the instrument and the test bench under the combined action of the effective components in the solution, to promote the detachment of the nucleic acid from the adsorbed surface, and then further combine the shear force generated during the wiping process to further contaminate the nucleic acid. Dissociated from the adsorption surface and entered the Space Cleaner™ reagent system to achieve the effect of removing surface nucleic acid contamination. Space Cleaner™ can non-enzymatically hydrolyze nucleic acid strands under normal temperature conditions to achieve the effect of decontamination. The main principle of Bottle B is Ag^+ plasma nanoparticles can rapidly bind to the mercapto ($-\text{SH}$) on protease in bacteria, interfering with bacterial DNA synthesis. At the same time, the binding of Ag^+ and protein makes the protease inactive, leading to bacterial death. Ag^+ has a high redox potential, so after the death of bacteria, it will be released freely from the bacteria and produce bactericidal and antimicrobial effects from cycle to cycle.

Operation Instructions

First spray the Space Cleaner™ Bottle A on the target surface three times, wipe it with a clean paper towel, without washing with water; Bottle B repeats the above operation. It can also soak the hard-to-clean small parts in the Space Cleaner™ Bottle A and Bottle B for 1 minute, then rinse them with clean water and dry them before installing them.

Step1: Spray three times on various surfaces, no need to wait, dry with clean paper at once.

Step2: Spray three times on various surfaces, no need to wait, dry with clean paper at once.

Precautions

Wear safety glasses and disposable gloves when handling reagents; This product can be used in glass, ceramics, plastics, rubber, steel and metals; Avoid using Space Cleaner™ on electronic devices, such as electric dispensers or balance, spray sterile paper towels with Space Cleaner™ and wipe them quickly.

Storage & Term of validity

Room temperature, avoid direct sunlight, valid for 12 months.

Product Specifications

Cat. No.	Product Name	产品规格
PCN01-0500	Space Cleaner™	Bottle A: 300 mL; Bottle B: 300 mL